

## EN71 -1 Mouth Actuated Toys Durability Tester



### Product introduction

Uses: To prevent mouth-actuated toys or mouthpieces from being mistakenly inhaled by children and causing suffocation, suitable for testing the safety performance of child mouth-actuated toys

Standard configuration: 1 mouth-actuated toy durability tester, 1 power cord, 1 trachea, 2 sets of test fixture sleeves, 2 fuses

Feature:

The instrument satisfies the air volume of at least 295 cubic centimeters +/- 10 cubic centimeters within 5 seconds, and can alternately blow and suck toys 10 times.

Equipped with a high-precision pressure gauge and a gas fitting suitable for different toys, the gas flow can be controlled by the valve.

## Technical Parameters

Test Speed	120+/-10/min
Testing times	10 times
Gas-pressure meter	Digital, Units, Kpa, The positive and negative pressure display
Dimensions	46.5x60x25cm
Weight	30kg
Standards	ISO 8124-1 16 CFR 1500 ASTM F963 4.8 EN-71
Accessories	Tube gripper 2pcs
Power Supply	220V/50Hz, 4A

Please review video of the 13.8 kPa Mouth Actuated Durability Tester operation as below link:

<https://www.youtube.com/watch?v=wJGBOELvoss&t=74s>